

## PRECISION AID.

### 3. CLAIMS 4 PRINTS.

1. Precision Aid key slot units are designed to funnel downward the tip of the eye drop bottle into the slot to be held securely over the center of the eye until the eye drop is delivered.

2. Precision Aid key slot unit are also designed to protrude outward in front of the brow bar of the glass frame 1/4" so that the tip of the eye drop bottle has a safety clearance above the eye of 1/2" and the slot holds the bottle steady until the drop is delivered without contamination.

3. Precision Aid key slot units are designed to be anchored into the most precise positioning instrument ever invented to service the eyes consistantly. The drops are delivered to the center of the eye directly without touching anything eliminating any chance of contamination and infection.

10045045-011602