

HOWARD KAPPLE
1927 FAIRMEADOWS DR.
BETTENDORF, IOWA. 52722.

PRECISION AID.

ABSTRACT.

This device is mounted on a very familiar instrument the eye glass frame. For positioning to service the eye it is precise every time you put them on. A key slot is mounted to the underside of the brow bar above each eye . These key slots are angled outward in front of the brow bar so that the tip of the eye drop bottle clears the eyes safely. The key slots are centered exactly in the center of the eye top to bottom and side to side. The key slots are made of copper. The shape of the key slots are triangular at the top to funnel down the eye drop bottle into the key slot where the bottle will rest until the drop is delivered. The bottom half of the glass frame was removed along with the glass. This was to provide room for one had to hold down the lower eye lid when delivering drops from the bottle with the other hand.

2025-04-04 10:45:04